

MC770B/MC770I

# TECHNO MULTI-MASSAGE CHAIR OWNER'S MANUAL




salon

**K**  
KEYTON

Model 770B/770I

## RECLINERS



### GENERAL / BACK RECLINER

#### Back Recliner.

##### Manual Control



To operate the Back Recliner independently, keep these keys pressed until the Recliner has reached the desired position.

Once the key is released, the Recliner will automatically stop.


#### General Recliner.

##### Automatic Control



To recline the Back and the Footrest simultaneously press these keys for a short time (less than one second).

Once pressed, the Chair will recline to its maximum. To stop the reclining movement at any specific position, press any of these keys again.



### FOOTREST RECLINER

##### Manual Control



To operate the Footrest Recliner independently, keep one of these keys pressed until the Recliner has reached the desired position.

Once the key is released, the Recliner will stop automatically.

##### Automatic Control

To recline the Footrest independently, press one of these keys for a short time (less than one second).

Once pressed, the Footrest will recline to its maximum. To stop the reclining movement at any specific position, press any of these keys again.

## IMPORTANT SAFETY PRECAUTIONS

Please read these precautions carefully before using the Massage Chair.

Warnings and Cautions described here below are essential for safe use of the Chair.

Be sure to follow the instructions carefully.

### **DANGER**

To reduce the risk of electric shock:

- 1 Always unplug this appliance from the electrical outlet immediately after using and before cleaning.

### **WARNING**

To reduce the risk of burns, fire, electric shock or injury to persons:

- 2 An appliance should never be left unattended when plugged in. Unplug from outlet when not in use.
- 3 This product should not be used by or near children. To avoid possible injury to children or pets, always return Chair to an upright position and fold down the Leg Massager carefully after using it. Make sure nothing can be caught in the Footrest mechanism. **SEE FIG. 1**
- 4 Close supervision is necessary when this appliance is used by or near handicapped or disabled people.
- 5 Use this appliance only for its intended use as described in this manual. Do not use attachments or devices not recommended by manufacturer.
- 6 Any attempt to modify the Chair by adding any device may result in fire, electric shock or injury to people.
- 7 Never operate this appliance if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. In this case return the appliance to a service center for examination and repair.
- 8 Do not carry this appliance by supply cord or use cord as a handle.
- 9 Keep the cord away from heated surfaces.
- 10 Never drop or insert any object into any operating part.
- 11 Do not use outdoors. Keep the Chair out of places exposed to high humidity (saunas, swimming pools) and do not splash water into any part of the Chair.
- 12 Do not operate where aerosol (spray) products are being used or where oxygen is being administered.
- 13 To disconnect, turn all controls to the off position and then remove plug from outlet.
- 14 Do not use extension cords.
- 15 Always keep children and pets away from the Chair during operation, especially while the reclining system is working. **SEE FIG. 5**
- 16 Before reclining the Chair Back or the Footrest, make sure that there are no obstructions. **SEE FIG. 2**
- 17 After each use, be sure to press the **POWER** key to turn off the Chair. As a further caution, unplug the power cord from the outlet.

### **IMPORTANT**

- Do not exceed 15 minutes on one part of the body.
- Do not exceed 30 minutes of use per day.
- Maximum user weight is 265 pounds (120Kg).

**SAVE THESE INSTRUCTIONS • FOR HOUSEHOLD USE ONLY**

## IMPORTANT SAFETY PRECAUTIONS

Please read these precautions carefully before using the Massage Chair.  
Warnings and Cautions described here below are essential for safe use of the Chair.  
Be sure to follow the instructions carefully.



## SAFETY PRECAUTIONS

The following individuals are requested to consult their physician prior to using this chair in order to prevent accidents or physical damages.

- 18** Make sure that there is no obstacle behind the Chair when reclining it. If the Back hits a wall or a column, it may malfunction. Be sure to have a wide space behind the Chair.
- 19** This Chair was intended for people who enjoy good health and wish to stay fit. Consult with your physician before using the Massage Chair if any of the following circumstances applies to you: malignant tumors, pregnancy, heart disease, or any other health condition that might affect your normal well-being. Do not use on any area that has received medical treatment without consulting a physician prior to use.
- 20** Never put any part of the body between the Massage Heads. Never put any part of the body between the Leg Massager Set Lever and the Armrest, or between the Armrest and the Chair Back. Do not place hands in mechanical portion of the Leg Massager located under the seat. Do not use on children.
- 21** Do not use without clothing. Thin clothing increases the effectiveness of the massage, but exposing the skin directly to the Massage Chair can irritate the skin.
- 22** Do not use for massaging the head, abdomen, elbows or knees.
- 23** Do not stretch the back muscles for more than 15 minutes at one time. Use only at gentle setting when first using the Chair, and do not use for more than 5 minutes. You can gradually increase duration and intensity as you adjust to the Massage Chair.
- 24** Do not massage any one point for longer than 2 minutes.
- 25** Excessive massaging can overstimulate the muscles and nerves and result in an adverse effect. To prevent excessive massaging, the Massage Chair comes with an automatic 15-minute timer. After that time, the Massage Heads will return to the restored position and stop.
- 26** While using the Chair, if you start feeling sick or if the massage seems painful, stop use immediately. If the pressure seems too strong, the movement of the Massage Heads can be stopped by pressing the **STOP** key.
- 27** When the Chair Back is in a reclined position, do not sit on it or sit with your legs resting on the Headrest. Do not sit on the Leg Massager, as the Chair may fall over.
- 28** Do not allow children to play on or around the Chair, especially during operation. **SEE FIG. 6**
- 29** Do not allow children to use the Chair and keep it out of their reach.
- 30** Do not use the Chair after meals.
- 31** This Massage Chair was developed for normal household use only.
- 32** Before operating the Salton® Massage Chair make sure there is no obstruction in front or behind it that would impede the reclining system to extend properly. **SEE FIG. 4**
- 33** If when you first use the Chair you feel the massage on your back is too strong, use a towel, blanket or small pillow to soften it.

**SAVE THESE INSTRUCTIONS • FOR HOUSEHOLD USE ONLY**

## IMPORTANT SAFETY PRECAUTIONS

Please read these precautions carefully before using the Massage Chair.

Warnings and Cautions described here below are essential for safe use of the Chair.

Be sure to follow the instructions carefully.

The following individuals are requested to consult their physician prior to using this chair in order to prevent accidents or physical damages.

## SAFETY PRECAUTIONS

- 34** This device should not be used over swollen areas or skin eruptions. Do not use if you have any kind of pain. In that case, consult a physician.
- 35** When excessive pressure is applied to the Massage Heads, they may stop moving for safety. When moving up the Heads from the lowered position, the body weight may stop the Head movement. In such case, slightly lift the body and let the Heads move up.
- 36** Keep the Chair out of places exposed to high humidity and do not splash water into any part of the Chair.
- 37** Hold the plug when removing it. Do not pull at the cord.
- 38** Do not damage, alter, forcibly bend, pull, twist or knot the power cord, and do not place the cord under heavy objects or a chair leg.
- 39** Avoid contact with the Chair if your hands or feet are wet.
- 40** Check that the voltage corresponds to the one indicated on the Chair (only alternating current).
- 41** To reduce the risk of electric shock, this appliance has a polarized plug (one blade is wider than the other). This plug will fit in a polarized outlet only one way. If the plug does not fit fully in the outlet, reverse the plug. If it still does not fit, contact a qualified electrician to install the proper outlet. Do not modify the plug in any way. **SEE FIG. 7**
- 42** Do not sit on the Controller.
- 43** Do not treat the covering roughly:
  - a) Keep sharp or pointed objects away from the upholstery of the Chair. Be careful not to drop in ashes or matches on the Chair.
  - b) Exposure to direct sunlight can cause fading or a color change of the covering.
- 44** If the Chair functions abnormally, turn off the power immediately and have the Chair checked by an authorized service center.
- 45** Do not attempt to repair the Chair yourself.
- 46** The Chair Back, where the Massager Heads are located, is a very delicate and important section, so do not hit it in any way.

SAVE THESE INSTRUCTIONS • FOR HOUSEHOLD USE ONLY

## IMPORTANT SAFETY PRECAUTIONS


Please read these precautions carefully before using the Massage Chair.  
Warnings and Cautions described here below are essential for safe use of the Chair.  
Be sure to follow the instructions carefully.



### CAUTION


To reduce the risk of electric shock or fire:

Figure 1



To avoid injuring any person, pet or especially children with the Footrest, always check below before lowering the Footrest.

Figure 2




To avoid injuring yourself or any person with the back of the Chair, especially children, check before reclining the Chair Back.

Figure 3



To avoid injury, do not let children sit on the Armrest.

Figure 4




Before operating the Massage Chair, make sure there is no obstruction in front or behind it.

Figure 5




Always keep children and pets away from the Chair during operation, especially while the reclining system is working.

Figure 6



Do not allow children to play on or around the Chair, especially during operation.

Figure 7



Polarized Plug

SAVE THESE INSTRUCTIONS • FOR HOUSEHOLD USE ONLY



## CONTROL PANEL




Figure 8




If the LEDs of the three upper rows are turned on when starting the Chair (see keys: 1 to 9), a security code has been programmed. Follow the instructions described in **Figure 9**.




If when you first use the Chair you feel the massage on your back is too strong, you can use a towel, a blanket or a small pillow to soften it, or request additional internal cushions from your distributor. Please see drawing below for instructions on how to place these internal cushions.

## INTERNAL CUSHIONS



For Higher Intensity



For Lower Intensity



### HOW TO TIE INTERNAL CUSHION



### NOTE:

Please untie the internal cushion before taking it off. Notice that there are four knots on both sides of the cushion. You must untie the eight knots. ■

After placing in an internal cushion, make sure that the eight knots are well tied up.



The lowest knot on both sides can be located inside a small opening in the upholstery. ■

## PERSONAL SECURITY CODE

To program the security code press four keys from the three upper rows of the Control Panel (numbered 1 to 9 according to Figure 8).

### USING THE CHAIR WHEN IT HAS A SECURITY CODE




- 1  Turn on the Chair by pressing **POWER**. (The LEDs of the keys 1-9 will blink)
- 2  Enter your personal security code.
- 3  The Chair works normally.

Figure 9

### PROGRAMMING THE SECURITY CODE





- 1  Before turning on the Chair, hold down **STOP** and then press **POWER**. (The LEDs of the keys 1-9 will blink)
- 2  The **LEGS** key will be turned on if the Chair does not have any memorized code. If it is not turned on, go to step 2 of figure 11.
- 3  Enter a security code.
- 4  The Chair works normally with a programmed security code.

Figure 10

### CHANGING THE SECURITY CODE







- 1  Before turning on the Chair, hold down **STOP** and then press **POWER**. (The LEDs of the keys 1-9 will blink)
- 2  The **LEGS** key is not turned on when there is a memorized security code. If the **LEGS** key is turned on, then the Chair does not have any programmed code.
- 3  Enter the security code using the keys 1-9 according to Figure 8.
- 4  The **LEGS** key will be turned on if the entered code is right.
- 5  Enter the new code.
- 6  The Chair works normally with a programmed security code.

Figure 11

### DELETING THE SECURITY CODE







- 1  Before turning on the Chair, hold down **STOP** and then press **POWER**. (The LEDs of the keys 1-9 will blink)
- 2  The **LEGS** key is not turned on when there is a memorized security code. If the **LEGS** key is turned on, then the Chair does not have any programmed code.
- 3  Enter the security code using the keys 1-9 according to Figure 8.
- 4  The **LEGS** key will be turned on if the entered code is right.
- 5  Press **STOP** to delete the security code.
- 6  The Chair works normally without a programmed security code.

Figure 12




### PLEASE REMEMBER YOUR SECURITY CODE


If you forget your personal security code, contact an authorized service center.





# MASSAGE PROGRAMMING / MEMORY A-B

## MEMORIZING MESSAGE SESSION "A"

- 

Before turning on the Chair, hold down **AUTO** and then press **POWER**.
- 


The LED of the **POWER** key will blink while the memorized program is being recorded.
- 


You can use the Chair normally. (The Chair will memorize all the functions selected.)
- 


Press the **STOP** key to stop recording the message session.


**Figure 13**


## MEMORIZING MESSAGE SESSION "B"

- 

Before turning on the Chair, hold down **AUTO** and then press **POWER**.
- 

The LED of the **POWER** key will blink while the memorized program is being recorded.
- 

Press **AUTO** to select memory B.
- 


You can use the Chair normally. (The Chair will memorize all the functions selected.)
- 


Press the **STOP** key to stop recording the message session.

Press the **STOP** key to stop recording the message session.


**Figure 14**

## STARTING THE PROGRAMS ALREADY MEMORIZED


- 

Turn on the Chair by pressing **POWER**.
- 

**PROGRAM A:**



**PROGRAM B:**

Press **AUTO** to start the memorized program A or hold down **AUTO** for over 3 seconds to start the memorized program B.
- 


The LED of the **AUTO** key will blink while the program lasts.

**Figure 15**

The Chair memorizes all the functions during 15 minutes or until the **STOP** key is pressed. Memory A can store up to 24 taps on the control panel, and memory B up to 25 taps.

The Chair comes with a default automatic program session memorized in programs A and B. If you want to restore this default automatic program after you have programmed them, select the program and press **STOP**.

## OPERATING THE CHAIR




### POWER ON/OFF

This key turns on or off the Chair power.  
The Chair should always be turned off whenever the massage programs are not in use.




If the LEDs of the three upper rows are turned on when starting the Chair (see keys 1 to 9), a security code has been programmed. Follow the instructions described in **Figure 9**.



### AUTOMATIC MEMORY PROGRAMS


This key initiates the automatic massage programs.  
See Figure 13 for instructions to memorize up to two automatic programs (A or B).



### MESSAGE STOP


This key stops all massage functions and brings the internal mechanism to its resting position.

## ADJUSTMENTS



### WIDTH ADJUSTMENT

By pressing this key repeatedly you can change the width of the massage surface.



### HEIGHT ADJUSTMENT

By pressing this key repeatedly you can adjust the reach of the massage according to your height.

## LOCATORS / LOCALIZED MESSAGE



### MESSAGE UP

By keeping this key pressed you can localize the massage on an upper area. By releasing this key, the vertical massage movement is localized on the desired point, and every two minutes it performs a partial massage.

You can initiate the partial massage with a short press (less than 1 second). In order to activate the complete vertical movement again, press any massage key (Figure 8: keys 2, 5, 8).




### MESSAGE DOWN

By keeping this key pressed you can localize the massage on a lower area. By releasing this key, the vertical massage movement is localized on the desired point, and every two minutes it performs a partial massage.


You can initiate the partial massage with a short press (less than 1 second). In order to activate the complete vertical movement again, press any massage key (Figure 8: keys 2, 5, 8).

## MASSAGES




### KNEADING/SENSITIVE KNEADING MESSAGE

By pressing this key repeatedly you will activate the different Kneading Massage programs.




### TAPPING MESSAGE

By pressing this key repeatedly you will activate the different Tapping Massage programs.



### ROLLING MESSAGE

By pressing this key the Rolling Massage system will be activated. This massage adapts itself to each natural curve of the back.



### LEG MESSAGE

By pressing this key repeatedly you will change the intensity of the massage. Once you have selected the intensity, you will have the possibility of initiating another Leg Massage program by pressing this key for 3 seconds. To change the program, press this key again shortly; to change the intensity, press the key for 3 seconds.

## RECLINERS

### GENERAL / BACK RECLINER

#### Back Recliner.

##### Manual Control



To operate the Back Recliner independently, keep these keys pressed until the Recliner has reached the desired position.

Once the key is released, the Recliner will automatically stop.

#### General Recliner.

##### Automatic Control



To recline the Back and the Footrest simultaneously press these keys for a short time (less than one second).

Once pressed, the Chair will recline to its maximum. To stop the reclining movement at any specific position, press any of these keys again.

### FOOTREST RECLINER

##### Manual Control



To operate the Footrest Recliner independently, keep one of these keys pressed until the Recliner has reached the desired position.

Once the key is released, the Recliner will stop automatically.

##### Automatic Control

To recline the Footrest independently, press one of these keys for a short time (less than one second).

Once pressed, the Footrest will recline to its maximum. To stop the reclining movement at any specific position, press any of these keys again.

1. The first part of the document is a list of the names of the persons who have been appointed to the various positions of the Board of Directors of the Corporation.

2. The second part of the document is a list of the names of the persons who have been appointed to the various positions of the Board of Directors of the Corporation.

3. The third part of the document is a list of the names of the persons who have been appointed to the various positions of the Board of Directors of the Corporation.

4. The fourth part of the document is a list of the names of the persons who have been appointed to the various positions of the Board of Directors of the Corporation.

5. The fifth part of the document is a list of the names of the persons who have been appointed to the various positions of the Board of Directors of the Corporation.

6. The sixth part of the document is a list of the names of the persons who have been appointed to the various positions of the Board of Directors of the Corporation.

7. The seventh part of the document is a list of the names of the persons who have been appointed to the various positions of the Board of Directors of the Corporation.

8. The eighth part of the document is a list of the names of the persons who have been appointed to the various positions of the Board of Directors of the Corporation.

9. The ninth part of the document is a list of the names of the persons who have been appointed to the various positions of the Board of Directors of the Corporation.

10. The tenth part of the document is a list of the names of the persons who have been appointed to the various positions of the Board of Directors of the Corporation.

11. The eleventh part of the document is a list of the names of the persons who have been appointed to the various positions of the Board of Directors of the Corporation.

12. The twelfth part of the document is a list of the names of the persons who have been appointed to the various positions of the Board of Directors of the Corporation.

13. The thirteenth part of the document is a list of the names of the persons who have been appointed to the various positions of the Board of Directors of the Corporation.