


Relaxor®



MOOD THERAPY LAMP

OWNER'S MANUAL &
OPERATING INSTRUCTIONS

Model ML8EGG

LIGHT THERAPY

IMPORTANT SAFEGUARDS

Warning — To Reduce the Risk of Fire, Burns, or Injury to Persons or Eye Damage:

- 1. NEVER TOUCH THE LAMP OR THE BASE DURING NORMAL USE. LIGHTED LAMP IS HOT!**
- 2. Disconnect Adapter from wall outlet and allow to cool for at least 25 minutes before touching the Lamp Shade or Base.*
- 3. Bulb gets HOT quickly! Be careful not to come into contact with the Lamp Shade or Base during use. Remove electric cord from wall to disconnect.*
- 4. Do not direct light at persons.*
- 5. Do not look directly at lighted Bulb.*
- 6. Keep Lamp away from curtains, draperies, and similar materials.*
- 7. This Lamp is designed for use with Adapter provided: 100-120V AC 60 hz. input and 12V AC output. Use of any adapter not recommended or provided by the manufacturer could result in injury to persons or property.*
- 8. Do not touch the Bulb at any time. Use a soft cloth. Oil from skin may damage the Bulb.*
- 9. Never allow water or other liquids to come into contact with the unit.*
- 10. Unplug unit from electrical outlet when not in use.*
- 11. For indoor use only.*

SAVE THESE INSTRUCTIONS FOR HOUSEHOLD USE ONLY

CAUTION


LIGHTED LAMP IS HOT!
Do not touch hot Bulb or Bulb Holder.

CONGRATULATIONS ON OWNING THE RELAXOR® MOOD THERAPY LAMP.

THE MOOD THERAPY LAMP OFFERS THE ULTIMATE IN
ENVIRONMENT-ENHANCING RELAXATION.

TO ASSEMBLE LAMP

1. Carefully unpack all of the Lamp components, and set up Lamp Base on a flat, stable surface.
2. Place the Shade on the Base.
3. Be sure that the Adapter is attached to the power cord.
4. Plug the Adapter into a standard 120V AC electrical outlet.
5. Turn the Switch on cord to light Mood Lamp.
6. To activate the Color Swirl, press Switch at Base of Lamp to ON, indicated by a "I".
7. When finished turn Switch on cord to deactivate Lamp.



USER MAINTENANCE INSTRUCTIONS

CLEANING

1. Unplug unit and allow it to cool.
2. Take a slightly damp cloth and wipe unit clean. Never allow water to come into contact with the Color Switch. Do not use any chemicals or abrasive cleaners.
3. With a dry cloth, wipe excess moisture off unit.
4. Do not attempt to clean Color Wheel or any other components inside the Base. Never immerse Base in water or other liquid.
5. Periodically untangle cord if it becomes twisted.

STORAGE

Unplug unit from wall and allow it to cool. Store unit in its box or in a clean, dry place.

LIMITED ONE YEAR WARRANTY

Warranty: This Relaxor® product is warranted by Salton, Inc. to be free from defects in materials or workmanship for a period of (1) year from the original purchase date. This product warranty covers only the original consumer purchaser of the product.

Warranty Coverage: This warranty is void if the product has been damaged by accident in shipment, unreasonable use, misuse, neglect, improper service, commercial use, repairs by unauthorized personnel, normal wear and tear, improper assembly, installation or maintenance abuse or other causes not arising out of defects in materials or workmanship. This warranty is effective only if the product is purchased and operated in the USA, and does not extend to any units which have been used in violation of written instructions furnished with the product or to units which have been altered or modified or to damage to products or parts thereof which have had the serial number removed, altered, defaced or rendered illegible.

Implied Warranties: ANY IMPLIED WARRANTIES WHICH THE PURCHASER MAY HAVE ARE LIMITED IN DURATION TO ONE (1) YEAR FROM THE DATE OF PURCHASE. Some states do not allow limitations on how long an implied warranty lasts, so the above limitation may not apply to you.

Warranty Performance: During the above one-year warranty period, a product with a defect will be either repaired or replaced with a reconditioned comparable model (at our option) when the product is returned to the Repair Center, or the purchase price refunded. The repaired or replacement product will be in warranty for the balance of the one-year warranty period and an additional one-month period. No charge will be made for such repair or replacement.

Service and Repair: Should the appliance malfunction, you should first call toll-free 1 (800) 233-9054 between the hours of 9:00 am and 5:00 pm Central Standard Time and ask for CONSUMER SERVICE stating that you are a consumer with a problem. Please refer to model number ML8EGG when you call.

In-Warranty Service (USA): For an appliance covered under the warranty period, no charge is made for service or postage. Call for return authorization 1 (800) 233-9054.

Out-of Warranty Service: A flat rate charge by model is made for out-of-warranty service. Include \$ 6.00 (U.S.) for return shipping and handling. We will notify you by mail of the amount of the charge for service and require you to pay in advance for the repair or replacement.

For Products Purchased in the USA, but Used in Canada: You may return the product insured, packaged with sufficient protection, and postage and insurance prepaid to the USA address listed below. Please note that all customs duty / brokerage fees, if any, must be paid by you and we will require you to pay the cost of customs duty / brokerage fees to us in advance of our performing any service.

Risk During Shipment: We cannot assume responsibility for loss or damage during incoming shipment. For your protection, carefully package the product for shipment and insure it with the carrier. Be sure to enclose the following items with your appliance: any accessories related to your problem, your full return address and daytime phone number, a note describing the problem you experienced, a copy of your sales receipt or other proof of purchase to determine warranty status. C.O.D. shipments cannot be accepted.

To return the appliance, ship to:

ATTN: Repair Center
Salton, Inc.
708 South Missouri Street
Macon, MO 63552

To contact us, please write to, call, or email:

Salton, Inc.
P.O. Box 1637
Columbia, MO 65205-1637
1-800-233-9054
E-mail: Salton@Saltonusa.com

Limitation of Remedies: No representative or person is authorized to assume for Salton, Inc. any other liability in connection with the sale of our products. There shall be no claims for defects or failure of performance or product failure under any theory of tort, contract or commercial law including, but not limited to negligence, gross negligence, strict liability, breach of warranty and breach of contract. Repair, replacement or refund shall be the sole remedy of the purchaser under this warranty, and **in no event shall Salton be liable for any incidental or consequential damages, losses or expenses.** Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

Legal Rights: This warranty gives you specific legal rights, and you may also have other rights which vary from state to state.

For more information on Salton products: visit our website: www.salton-maxim.com

©2001 Salton, Inc. All rights reserved.

Made in China

P/N 60808

8/01